

# MitraClip

## Percutaneous Mitral Valve Repair

### Introduction to the MitraClip Therapy

# Agenda

Mitral Regurgitation

Concept & Device Overview

MitraClip Procedure

Clinical experience

Patient profiles

# Mitral Regurgitation

# About Mitral Regurgitation

## Mitral Regurgitation

Mitral regurgitation (MR) results from a lack of leaflet coaptation between the two leaflets of the mitral valve—the valve found between the left atrium and left ventricle of the heart. In a normally functioning mitral valve, blood flows in a single direction between the left atrium and left ventricle. MR is characterized by systolic retrograde flow from the left ventricle into the left atrium. Over time, MR may lead to heart failure.

## Causes

- Degenerative MR (also known as primary or organic MR) is usually due to an anatomic abnormality of the mitral valve itself, including the leaflets, and/or the subvalvular apparatus, such as the chordae or papillary muscles.
- Functional MR (also known as secondary MR) is the result of left ventricular dilation, which can be secondary to ischemic heart disease. Left ventricular dysfunction leads to annular dilation and incomplete coaptation of the mitral valve resulting in MR.

## Consequences

MR places an extra burden on the heart and lungs. In some cases, patients may never develop symptoms. In other cases, patients may develop an enlarged left ventricle as the heart works harder to maintain cardiac output. Patients with chronic MR may develop symptoms of heart failure, such as fatigue or inability to exercise, decrease in appetite, dry, hacking cough (often worse when lying down), shortness of breath especially at night, fainting, weight gain from fluid retention, or accumulation of fluid in feet, ankles, and lungs (edema).

Patients can also develop rhythm abnormalities that may require medical or interventional therapy. Overall patients with symptomatic severe MR have a worse quality of life, decreased functional capacity and increased rate of mortality.



Degenerative MR — Prolapse



Degenerative MR — Flail





Functional MR

# Prevalence and Unmet Need

Epidemiological data show that MR in the range of moderate-to-severe to severe is the most frequent valve disease.<sup>1,2</sup> In fact, MR is the second most common form of valvular heart disease needing surgery in Europe.<sup>3</sup>

Although moderate-to-severe and severe MR are common, they continue to remain undertreated. In 2009, only approximately 2% of the estimated 1.7 million patients with MR  $\geq$  3+ were treated with surgical intervention.<sup>1,7-10</sup> Reasons for denying surgery include impaired left ventricular ejection fraction (LVEF), multiple comorbidities, and advanced age, all of which are determinants for high operative risk.<sup>3,6</sup>



1 Heart Disease and Stroke Statistics 2010 Update: A Report from the American Heart Association. Circulation 2010;121:e46-e215.

2 Nkomo, VT et al. Lancet. 2006;368:1005-1011.

3 Jung, B et al. Eur Heart J. 2003;24:1231-1243

4 Mirabel, M et al. Eur Heart J. 2007;28:1358-1365.

5 Taramasso, M et al. Cardiol Res Pract. 2010;25:[Epub ahead of print].

7 Rankin, J et al. Determinants of Operative Mortality in Valvular Heart Surgery. J Thorac Cardiovasc Surgery, March 2006.

8 U.S. Census Bureau, Statistical Abstract of the U.S.

9 Patel et al. Mitral Regurgitation in Patients with Advanced Systolic Heart Failure, J of Cardiac Failure, 2004.

10 Gammie, J et al. Trends in Mitral Valve Surgery in the United States: Results from the STS Adult Cardiac Database, Annals of Thoracic Surgery 2009.

11 Data on file Abbott Vascular, March 2011.

# Treatment Options

**Medical:** There are no medications that are indicated to treat MR, but there are medications used to manage patient symptoms.

**Surgical:** For symptomatic patients diagnosed with moderate–severe or severe MR, surgery is generally recommended to repair or replace the mitral valve. Mitral valve repair or replacement typically involves open-heart surgery while on cardiopulmonary bypass. Patients recovering from mitral valve surgery may take several months to regain normal physical function and activity.


**Percutaneous mitral valve repair:** The MitraClip procedure is a minimally invasive catheter-based therapy. This new treatment increases the options for selected patients with MR. It has been shown to reduce MR, reverse left ventricular remodeling, improve NYHA functional class, and improve quality of life.<sup>5</sup>




5. Feldman T et al. EVEREST II RCT, TCT 2012, NPL03998-A

# Concept

# Concept: Percutaneous Mitral Valve Repair



- Double-orifice suture technique developed by Prof. Ottavio Alfieri
- First published results in 1998 illustrated proven benefit
- Suggested procedure best suited for minimally invasive approach



- Dr. Fred St. Goar, interventional cardiologist had patient successfully treated with edge-to-edge surgery
- Conceived several ideas for percutaneous valve repair
- Founded Evalve 1999 to develop devices to treat valvular disease

# Device Overview

# System

The MitraClip System performs **percutaneous mitral valve repair** by creating a vertical line of coaptation, forming a double-orifice valve.

- Beating heart procedure—no cardiopulmonary bypass
- Allows for real-time positioning and repositioning to optimize MR reduction
- Designed to preserve surgical options
- Femoral venous access
- Limited hospital length of stay compared to that after surgery<sup>12</sup>



The Clip Delivery System and Steerable Guide Catheter are designed to fit co-axially to accurately position and reposition multiple implants with the use of one guide.

12. Data on file Abbott Vascular. November 2011

# System Components

## Steerable Guide Catheter

- 24 French steerable catheter
- Percutaneous venous access
- Requires successful transeptal puncture


## Clip Delivery System

Contains the implant attached to a highly maneuverable delivery catheter with all controls at the proximal end.

## MitraClip Device (Implant)

- Cobalt chromium construction
- Polyester cover designed to promote tissue growth
- Magnetic resonance conditional to 3 Tesla\*

\* Static magnetic field up to 3 Tesla; maximum spatial gradient in static field of 2500 gauss/cm or less; maximum whole-body averaged specific absorption rate (SAR) of 3.0 W/kg for 15 minutes of scanning.




# MitraClip Procedure


# MitraClip Procedure Animation




# Cut Chord Model




# Procedural Overview




Venous access under general anesthesia



Clip is optimally positioned on MR jet




Creation of double-orifice valve



Clip implanted after validation of safe positioning

# Procedural Overview


## Patient and System Preparation




- The following considerations should be accounted for:
  - TEE probe will be in place for an extended period of time
  - Intubation under general anesthesia
  - 24 French sheath in femoral vein
  - Bladder/urinary catheter in place
  - Heparinization during procedure to ACT > 250
- System is prepared by removing all the air in the lumens of the Clip Delivery System and Steerable Guide Catheter
- System is functionally tested prior to use

# Procedural Overview

## Transseptal Crossing and Guide Insertion




A transseptal procedure is performed to gain access from the right atrium to the left atrium. The Steerable Guide Catheter (Guide) and Dilator are then carefully advanced into the left atrium over a wire. Once the Guide is in place and secured, the wire and Dilator are removed leaving the Guide in the left atrium.





# Procedural Overview

## Clip Delivery System Insertion and Steering in the Left Atrium




To introduce the Clip, the Clip Delivery System (CDS) is advanced through the Guide into the left atrium. A series of steering maneuvers and manipulations with the Guide and CDS are required to align the Clip perpendicular to the mitral valve plane, and the Clip Arms perpendicular to the line of coaptation. These maneuvers are done under echocardiographic and fluoroscopic guidance.




# Procedural Overview

## Advancing into Left Ventricle and Leaflet Grasping



After the Clip is aligned over the regurgitant jet in the left atrium, the System is then advanced into the left ventricle to begin the grasping procedure. Leaflet grasping is done by slowly retracting the System back towards the left atrium to allow the leaflets to come to rest on the Clip Arms and then dropping the Grippers.




Clip Arms closed to 120°

Clip Arms closed to 20°




# Procedural Overview

## Leaflet Insertion Assessment and Hemodynamic Measurements




Prior to Clip closure and deployment, a leaflet insertion and hemodynamic assessment must be performed. The leaflet insertion assessment ensures both leaflets are fully inserted and secure into the Clip. In addition, the MR reduction and pressure gradients are assessed to ensure regurgitation reduction without stenosis.



# Procedural Overview

## Deployment and System Removal




Once the assessments are positive, the Clip can be fully closed and deployed in a multistep process. The physician may also decide to place a second Clip to optimize MR reduction. The System is removed by releasing deflections on the catheter and slowly removing from the patient. Groin management and continued medical therapy are recommended per the institution's guidelines.



# Post-procedure Considerations

## Post-Procedure Recovery Instructions

- Patient might have to be intubated, procedure performed under general anesthesia
- Patient may have Femoral Arterial and/or Venous access
- Patient will have had 24 French sheath in Femoral Vein during procedure
- Patient will have foley catheter in place
- Patient will have had TEE probe in place for extended period of time
- Antibiotic Therapy
  - Administer prophylactic antibiotics per institutional guidelines for implanted devices
- Groin Access
  - Per institutional guidelines and similar to other catheterization procedures
- Anticoagulation Therapy
  - Short-term anticoagulation therapy may be necessary after cardiac valve repair with the MitraClip device. Prescribe anticoagulation and other medical therapies per institutional guidelines.



# Clinical Experience

# Worldwide Clinical Experience

- Over 13 000 patients have been treated with the MitraClip Therapy worldwide.<sup>1</sup>
  - 75% are considered high risk\* for mitral valve surgery
  - 67% have functional mitral regurgitation (MR)
- The use of the MitraClip is supported by a rigorous clinical trial program.<sup>1</sup>
  - 50% are considered high risk\* for mitral valve surgery




1. Data as of March 2014

\*Determination of high surgical risk based on: logistic EuroSCORE  $\geq$  20%, or STS calculated mortality  $\geq$  12%, or pre-specified high surgical risk co-morbidities specified in EVEREST II High Risk Study protocol.

# MitraClip Therapy

## Current Global Adoption




<b>Treating Centers</b>	328
<b>Patients</b> (clinical and commercial)	12,431
<b>Patients</b> <sup>1</sup> (commercial)	11,089
<b>Implant Rate</b> <sup>2</sup>	96%
<b>Acute MR reduction</b> <sup>2</sup>	99% of implants
<b>Functional MR</b> <sup>2</sup>	67%
<b>Degenerative MR</b> <sup>2</sup>	21%
<b>Mixed</b> <sup>2</sup>	12%

■ Commercial      ■ Special Access

1. First-time procedures only. Includes commercial patients, ACCESS I and ACCESS II patients.  
 2. Successful implants only.

# Growing body of clinical evidence




\*large cohort is defined as >50 patients

# MitraClip Therapy

## Broad Spectrum of Experience

### EVEREST II


(Randomized Controlled Trial)



- 178 patients
- Device time – 156 minutes
- Implant rate – 89%

### EVEREST


(High Risk Cohort<sup>^</sup>)



- 211 patients
- Device time – 128 minutes
- Implant rate – 95%

### ACCESS EU

(Europe)



- 567 patients
- Device time – 118 minutes
- Implant rate – 99.6%

■ = DMR    ■ = FMR

<sup>^</sup> Enrolled by February 28, 2010  
Data on file Abbott Vascular, April 12, 2011  
Schillinger, W. ESC 2012, ACCESS 1-year results


Prospective, randomized, parallel-controlled, multi-center clinical evaluation of the MitraClip device plus optimal standard of care therapy (Device group) compared to optimal standard of care therapy alone (Control group).

- Randomized 1:1 between MitraClip + standard of care therapy or standard of care therapy alone
- Approximately 800 patients to be enrolled in up to 75 sites across Europe
- Trial device: MitraClip system used within the IFU specifications
- 24mo follow up for all patients

# Patient profile

# Multidisciplinary Team

A multidisciplinary approach and collaboration across specialties is critical to MitraClip Therapy success.



# MitraClip in Specific Patient Populations

## Patient groups in which significant clinical benefits have been reported:

- Degenerative MR, declined for surgery<sup>1</sup>
- Severe LV dysfunction refractory to medical therapy<sup>2</sup>
- Severe Heart Failure, despite optimal medical therapy<sup>3</sup>
- CRT non-responders<sup>4</sup>
- Bivalvular Disease: Severe Aortic Stenosis and Mitral Regurgitation<sup>5</sup>

## The following parameters should be taken into consideration by the Heart Team<sup>6</sup>:

- Moderate to severe or severe MR (Functional or Degenerative)
- Echocardiographic criteria for eligibility
- Level of surgical risk
- Greater than one year life expectancy


1.Reichenspurner, H. et al. Clinical Outcomes through 12 months in patients with Degenerative Mitral Regurgitation treated with the MitraClip device in the ACCESS-Europe Phase I trial. Eur J Cardiothoracic Surgery. 2013; 44:e 280-288. 2. Franzen O, Baldus S, Rudolph V, et al. Acute outcomes of MitraClip therapy for mitral regurgitation in high-surgical-risk patients: Emphasis on adverse valve morphology and severe left ventricular dysfunction. Eur Heart J. 2010; 31:1373-1381. 3. Franzen et al. MitraClip Therapy In Patients With End-Stage Systolic Heart Failure. Eur J Heart Failure. 2011; 13: 569-576. 4. Auricchio et al. Correction of Mitral Regurgitation in Nonresponders To Cardiac Resynchronization Therapy By MitraClip Improves Symptoms And Promotes Reverse Remodeling. JACC 2011; 58: 2183-2189. 5. Rudolph V, Schirmer J, Franzen O, Schlüter M, Seiffert M, Treede H, Reichenspurner H, Blankenberg S, Baldus S. Bivalvular transcatheter treatment of high-surgical-risk patients with coexisting severe aortic stenosis and significant mitral regurgitation. Int J Cardiol. 2013; 167(3):716-20. 6. ESC/EACTS 2012 Guidelines on the management of valvular heart disease. Eur Heart J (2012) 33, 2451–2496.  
The data is not from prospective studies and study results should be interpreted with caution





# MitraClip “beating porcine model” (for demonstration only)

# Porcine beating heart model / MitraClip procedure



**Abbott Vascular International BVBA**

Park Lane, Culliganlaan 2B, B-1831 Diegem, Belgium, Tel: 32.2.714.14.11

EVEREST I, EVEREST II, REALISM, ACCESS-EU, COAPT & RESHAPE-HF are Abbott Vascular Sponsored Studies.

Product is subject to prior training requirement as per the Instruction for Use. This product is intended for use by or under the direction of a physician. Prior to use, it is important to read the package insert thoroughly for instructions for use, warnings and potential complications associated with the use of this device. Information contained herein is for distribution for Europe, Middle East and Africa ONLY. Please check with the regulatory status of the device before distribution in areas where CE marking is not the regulation in force.

All drawings are artist's representations only and should not be considered as an engineering drawing or photograph. Photo(s) on file at Abbott Vascular.

The following needs to be considered by French healthcare professionals only.

Clip de réparation mitrale MitraClip et accessoires. Dispositifs médicaux de classe III et I, organisme notifié Dekra. Fabriqué par Evalve Inc, mandataire européen MedPass. Se référer aux informations de la notice d'instructions qui décrivent les informations de bon usage du dispositif. Veuillez lire attentivement les instructions figurant dans la notice. Non pris en charge par les organismes d'assurance maladie.

MitraClip is a trademark of the Abbott Group of Companies.

For more information, visit our website at [abbottvascular.com](http://abbottvascular.com)

© 2014 Abbott. All rights reserved. 9-EH-4-4331-01 04-2014

